Technical Sheet

GL SPECTIS 1.x series

Precise light measurement technology in a handy size.

GL **SPECTIS 1.x** is a high quality, easy to operate measuring device that gives you all you need for reliable light measurement. It is the perfect instrument for the measurement of LEDs, as well as for the final assessment of lamps or testing of complete lighting installations.

Features:

- High sensitivity and precise calibration
- Low noise and stable measurements
- One device can work with multiple optical probes
- Ready to work when connected to PC
- Small size and low energy consumption
- Powered via USB connection
- Available in options for different spectral and sensitivity ranges



Natural light, LEDs, halogen light, el	c.	
LED MEASUREMENT		
llluminance (lux)*	10 lx 100 000 lx 10 lx 100 000 lx - 5 lx 50 000 lx	Spectis 1.0 Spectis 1.1 Spectis 1.2 Spectis 1.3
Luminance [cd/m²]	Available with optional GL Opti Probe	
Luminous flux [lm]	Available with optional GL Opti Sphere	
Luminous intensity [cd]	Calculated in GL Spectrosoft	
Illuminance class	Class B – DIN 5032-7 Class AA – JIS C 1609-1:2006	
Tolerance – cosine response (f2')	< 3 % (1.9 %)	
Spectral range**	340 – 780 nm (UVa – VIS) 340 – 780 nm (UVa – VIS) 640 – 1050 nm (VIS – NIR) 340 – 750 nm (UVa – VIS)	Spectis 1.0 Spectis 1.1 Spectis 1.2 Spectis 1.3

GL OPTIC Polska Sp. z o.o. Sp.k

ul. Poznańska 70, 62-040 Puszczykowo, Poland Phone: +48 61 819 40 03 | E-mail: office@gloptic.com www.gloptic.com



Technical Sheet

GL SPECTIS 1.x series

Precise light measurement technology in a handy size.

Sensor		CMOS image sensor
Number of pix	els	256
	ition / datapoint interval	~ 1.7 nm
Wavelength re		0.5 nm
Integration tim		5 ms – 10 s
A/D converter		16 bit
Signal-to-noise	e ratio	1000:1
Stray light		2*10 E-3
Optical resolut	ion / FWHM	10nm
	curacy ***/****	5 % within range 340 – 500 nm 4 % within range 500 – 1050nm
Flicker comper	nsation	√
Temperature s and dark curre	ensor ent compensation	\checkmark
Uncertainty of	color coordinates ***	0.0015
Automatic acce	essory detection	\checkmark
GENERAL PE	ROPERTIES	
Power supply	via USB connector	< 640 mA
Operating tem	perature	5 – 35 °C
Dimensions [H	× W × D]	62 mm x 115 mm x 28,3 mm (with standard diffuser)
Weight		125 g
Tripod adapter		\checkmark
INTERFACE	& MEMORY	
USB		USB 2.0
Measurement	result storage	by GL Spectrosoft
Trigger		Available with GL SPECTIS 1.1 Open collector, minijack 3.5mm, 3-pin
Data format		XML
SOFTWARE		
Software		Optional GL SPECTROSOFT Basic / Pro / Lab
ORDERING I	NFORMATION	
Case		\checkmark
USB cable		✓
Part number:	Spectis 1.0 VIS Spectis 1.1 VIS Spectis 1.2 VIS-NIR Spectis 1.3 LS	GLX10 no. 67827 GLX11 no. 106294 GLX12 no. 106302 GLX13 no. 202031

Dynamic range is spectrum related and should be calculated separately for any light source. Estimated dynamic range for typical 4000 K white LED.

Range estimated for optical system made to default specification. Alterations of that are often possible. Please consult technical support if you are looking for specific parameters. ** Spectral range of the sensor. Actual spectral range of system may be reduced due to limitations of used optical accessory.

*** Absolute measurement uncertainty immediately after calibration. The expanded uncertainty corresponds to a coverage probability of 95 % and the coverage factor k = 2. Parameters valid in laboratory conditions 25deg C, relative humidity 45%.

**** Applies only within the spectral range of the given model.

Note: Instrument, firmware and software specification are subject to change without prior notice. All information included in GL OPTIC datasheets and product information available in any form are carefully prepared and included information believed to be true. Please note that discrepancies may occur due to text and/or other errors or changes in the available technology. We advise to contact GL Optic before the use of the product to obtain the latest product specification.

GL OPTIC Polska Sp. z o.o. Sp.k



ul. Poznańska 70, 62-040 Puszczykowo, Poland Phone: +48 61 819 40 03 | E-mail: office@gloptic.com www.gloptic.com